

Promote

OSTEOPROTM DBM100

Part Numbers	Description	Size
PRO-01-05	OsteoPro TM 100% DBM	5cc
PRO-01-10	OsteoPro TM 100% DBM	10cc

Promote

OSTEOPROTM DBM100



AllianceSpineTM
Inspiring Solutions *Through Innovation*

14206 Northbrook Dr
San Antonio, TX 78232

P:(210) 314-2525 • F:(210) 314-2524
E: info@alliance-spine.com • www.alliance-spine.com
MKT-704 Rev. A




AllianceSpineTM
Inspiring Solutions *Through Innovation*

Promote

OSTEOPro™ DBM100

OsteoPro™ DBM100 is 100% dehydrated DBM, processed using proprietary demineralization process. It may be hydrated with saline, blood, Bone Marrow Aspirate (BMA), Platelet Rich Plasma (PRP), or other cellular components in accordance with a physician's well-informed medical judgment. **OsteoPro™ DBM100** does not contain any extrinsic carriers and is entirely derived from 100% human allograft bone.



Uses and Benefits of OsteoPro™ DBM100

- Bone filler for the treatment of surgically created or traumatic skeletal defects.
- 100% Bone – NO ADDITIONAL CARRIER
- Osteoinductive – each lot is assayed
- Osteoconductive
- Graft will not “wash away” during irrigation

OsteoPro™ DBM100 Advantages

- Stores dry at room temperature
- Suitable for immediate use off the shelf
- Sterilized to SAL 10⁻⁶ in accordance with ISO 11137
- Five year shelf life

Clinical Applications

- Spine /Neurologic Surgeries
- Orthopedic Surgery
- Foot & Ankle
- Reconstruction
- Trauma
- Dental



Tissue Safety

- Strict quality assurance and quality control procedures to ensure patient safety.
- Collection of donor tissue is under extensive medical review for medical /social history to determine eligibility
- Processed in accordance with the safety guideline provided by the U. S. the Food and Drug Administration (FDA) – Human Cellular and Tissue-based Products (HCT/P) (21 CFR Part 1271)
- Processed in accordance with the standards from the American Association of Tissue Banks (AATB)